

## AN ASSAY METHOD FOR SIMULTANEOUS ESTIMATION OF AMLODIPINE AND BENAZEPRIL IN ITS BULK AND PHARMACEUTICAL DOSAGE FORM BY RP-HPLC METHOD

Fatima Mazher\*<sup>1</sup> and Dr. R. Vani<sup>1</sup>

Department of Pharmaceutical Analysis, Shadan Women's College of Pharmacy, Hyderabad.

\*Corresponding Author: Fatima Mazher

Department of Pharmaceutical Analysis, Shadan Women's College of Pharmacy, Hyderabad.

Article Received on 20/09/2019

Article Revised on 10/10/2019

Article Accepted on 30/10/2019

### ABSTRACT

A short selective, precise, accurate and sensitive for the estimation of Benazepril and Amlodipine was done by RP-HPLC. The assay of Benazepril and Amlodipine was performed with tablets and the % assay was found to be 100.08 and 100.66 which shows that the method is useful for routine analysis. The linearity of Benazepril and Amlodipine was found to be linear with a correlation coefficient of 0.999 and 0.999, which shows that the method is capable of producing good sensitivity. The acceptance criteria of precision is RSD should be not more than 2.0% and the method show precision 0.8 and 0.9 for Benazepril and Amlodipine which shows that the method is precise. The acceptance criteria of intermediate precision is RSD should be not more than 2.0% and the method show precision 0.5 and 0.7 for Benazepril and Amlodipine which shows that the method is repeatable when performed in different days also.

**KEYWORDS:** Benazepril and Amlodipine, Validation, stability indicating method, degradation products.

### INTRODUCTION

Benazepril, 2-[(3S)-3-[[[(2S)-1-ethoxy-1-oxo-4-phenylbutan-2-yl]amino]-2-oxo-2,3,4,5-tetrahydro-1H-1-benzazepin-1-yl]acetic acid brand name Lotensin, is a medication used to treat high blood pressure (hypertension), congestive heart failure, and chronic renal failure. Upon cleavage of its ester group by the liver, benazepril is converted into its active form benazeprilat, a non-sulphydryl angiotensin-converting enzyme (ACE) inhibitor.




Fig. 1: Benazepril, 2-[(3S)-3-[[[(2S)-1-ethoxy-1-oxo-4-phenylbutan-2-yl] amino]-2-oxo-2,3,4,5-tetrahydro-1H-1-benzazepin-1-yl]acetic acid

Amlodipine, 3-ethyl 5-methyl 2-[(2-aminoethoxy) methyl]-4-(2-chlorophenyl)-6-methyl-1,4-dihydropyridine -3,5-dicarboxylate initially approved by the FDA in 1987, is a popular antihypertensive drug belonging to the group of drugs called *dihydropyridine calcium channel blockers*. Due to their selectivity for the peripheral blood vessels, dihydropyridine calcium channel blockers are associated with a lower incidence of myocardial depression and cardiac conduction abnormalities than other calcium channel blockers.




Fig. 2 Amlodipine, 3-ethyl 5-methyl 2-[(2-aminoethoxy) methyl]-4-(2-chlorophenyl)-6-methyl-1,4-dihydropyridine-3,5-dicarboxylate

**Experimental****Optimized chromatographic conditions**

Instrument used : Waters HPLC with auto sampler and UV detector.  
 Temperature : Ambient  
 Column : Inertsil ODSC18 (4.6×150mm 5µm)  
 Buffer : Phosphate buffer pH 4.8  
 pH : 3.0  
 Mobile phase : Buffer: Acetonitrile (35:65)  
 Flow rate : 0.8 ml per min  
 Wavelength : 242 nm  
 Injection volume : 20 µl  
 Run time : 10 min.

**Preparation of buffer and mobile phase****Preparation of phosphate buffer pH 4.8**

Weigh accurately 6.8gms of potassium di hydrogen orthophosphate dissolved in 1000 ml of HPLC water Ph was adjusted up to 3.0. Final solution was filtered through 0.44 µm Membrane filter and sonicate it for 10 mins.

**Preparation of mobile phase**

Accurately measured 350 ml (35%) of above buffer and 650 ml of Acetonitrile HPLC (65%) were mixed and degassed in an ultrasonic water bath for 10 minutes and then filtered through 0.45 µ filter under vacuum filtration.

**Diluent Preparation**

The Mobile phase was used as the diluent.

**Preparation of the benazepril & amlodipine standard & sample solution****Standard Solution Preparation**

Accurately weigh and transfer 40 mg of Benazepril and 20 mg of Amlodipine working standard into a 100 ml clean dry volumetric flask add about 7 mL of Diluent and sonicate to dissolve it completely and make volume up to the mark with the same solvent. (Stock solution).

Further pipette 1.5 ml of the above stock solutions into a 10ml volumetric flask and dilute up to the mark with diluent.

**Sample Solution Preparation**

Accurately weigh 10 tablets crush in mortar and pestle and transfer equivalent to 40 mg of Benazepril and 20 mg Amlodipine sample into a 100 mL clean dry volumetric flask add about 7 mL of Diluent and sonicate it up to 15 mins to dissolve it completely and make volume up to the mark with the same solvent. Then it is Filtered through 0.45 micron Injection filter. (Stock solution).

Further pipette 1.5ml of Benazepril and Amlodipine from the above stock solution into a 10ml volumetric flask and dilute up to the mark with diluent.

**Procedure**

Inject 20 µl of the standard, sample into the chromatographic system and measure the areas for Benazepril and Amlodipine peaks and calculate the % Assay by using the formulae.


**RESULTS AND DISCUSSION**

Fig. 3: Standard Chromatogram of Benazepril and Amlodipine.

Table 1: Results of system suitability parameters.

| S. No. | Name       | RT(min) | Area (µVsec) | Height (µV) | USP resolution | USP tailing | USP plate count |
|--------|------------|---------|--------------|-------------|----------------|-------------|-----------------|
| 1      | Benazepril | 2.227   | 121623       | 75723       |                | 1.54        | 3765.98         |
| 2      | Amlodipine | 5.819   | 42389        | 23721       | 7.98           | 1.38        | 5742.70         |

Table 2: Results of Linearity of Benazepril and Amlodipine.

| S. No. | Benazepril                         |        | Amlodipine                         |       |
|--------|------------------------------------|--------|------------------------------------|-------|
|        | Concentration ( $\mu\text{g/ml}$ ) | Area   | Concentration ( $\mu\text{g/ml}$ ) | Area  |
| 1      | 20                                 | 42486  | 10                                 | 14430 |
| 2      | 40                                 | 82149  | 20                                 | 29766 |
| 3      | 60                                 | 124321 | 30                                 | 43753 |
| 4      | 80                                 | 161418 | 40                                 | 57246 |
| 5      | 100                                | 201924 | 50                                 | 72471 |




Figure 4: Calibration graph for Benazepril.




Figure 5: Calibration graph for Amlodipine.

Table 3: Results of Precision for Benazepril.

| Injection                 | Area     |
|---------------------------|----------|
| Injection-1               | 123578   |
| Injection-2               | 122050   |
| Injection-3               | 123931   |
| Injection-4               | 121272   |
| Injection-5               | 123137   |
| Injection-6               | 122702   |
| <b>Average</b>            | 122778.3 |
| <b>Standard Deviation</b> | 990.2    |
| <b>%RSD</b>               | 0.8      |

**Table 4: Results of Precision for Amlodipine.**

| Injection                 | Area    |
|---------------------------|---------|
| Injection-1               | 43604   |
| Injection-2               | 43572   |
| Injection-3               | 43902   |
| Injection-4               | 43587   |
| Injection-5               | 44286   |
| Injection-6               | 44447   |
| <b>Average</b>            | 43899.7 |
| <b>Standard Deviation</b> | 385.1   |
| <b>%RSD</b>               | 0.9     |

**Table 5: Accuracy (recovery) data for Benazepril.**

| %Concentration<br>(at specification Level) | Area     | Amount Added<br>(mg) | Amount Found<br>(mg) | % Recovery | Mean<br>Recovery |
|--|----------|----------------------|----------------------|------------|------------------|
| 50%  | 60984.7  | 20                   | 20.04                | 100.21     | 100.42           |
| 100%                                       | 122049.7 | 40                   | 40.11                | 100.28     |                  |
| 150%                                       | 184002.3 | 60                   | 60.47                | 100.78     |                  |

\*Average of three determinations

**Table 6: Accuracy (recovery) data for Amlodipine.**


| %Concentration<br>(at specification Level) | Area    | Amount Added<br>(mg) | Amount Found<br>(mg) | % Recovery | Mean<br>Recovery |
|--|---------|----------------------|----------------------|------------|------------------|
| 50%  | 21510   | 10                   | 10.07                | 100.69     | 100.64           |
| 100%                                       | 43071.7 | 20                   | 20.16                | 100.81     |                  |
| 150%                                       | 64352.7 | 30                   | 30.12                | 100.41     |                  |

**Table 7: Results of LOD.**

| Drug name  | Baseline noise ( $\mu$ V) | Signal obtained ( $\mu$ V) | S/N ratio |
|------------|---------------------------|----------------------------|-----------|
| Benazepril | 54                        | 164                        | 3.04      |
| Amlodipine | 54                        | 162                        | 3.00      |

**Table 8: Results of LOQ.**

| Drug name  | Baseline noise ( $\mu$ V) | Signal obtained ( $\mu$ V) | S/N ratio |
|------------|---------------------------|----------------------------|-----------|
| Benazepril | 54                        | 541                        | 10.02     |
| Amlodipine | 54                        | 543                        | 10.06     |

**Degradation Studies****Figure 6: Chromatogram showing Acid degradation.**




Figure 7: Chromatogram showing Base degradation.




Figure 8: Chromatogram showing Peroxide degradation.




Figure 9: Chromatogram showing Thermal degradation.




Figure 10: Chromatogram showing Photo degradation.

Table 9: Results for Stability of Benazepril and Amlodipine.

| Sample Name | Benazepril |            | Amlodipine |            |
|-------------|------------|------------|------------|------------|
|             | Area       | % Degraded | Area       | % Degraded |
| Standard    | 121403     |            | 42673      |            |
| Acid        | 119736     | 1.37       | 41662      | 2.37       |
| Base        | 117533     | 3.19       | 41776      | 2.10       |
| Peroxide    | 113863     | 6.21       | 40623      | 4.80       |
| Thermal     | 116563     | 3.99       | 40572      | 4.92       |
| Photo       | 115552     | 4.82       | 39772      | 6.80       |

Table 10: Results of Assay for Benazepril and Amlodipine.

|            | Label Claim (mg) | % Assay |
|------------|------------------|---------|
| Benazepril | 20               | 100.08  |
| Amlodipine | 10               | 100.66  |

## SUMMARY AND CONCLUSION

The estimation of Benazepril and Amlodipine was done by RP-HPLC.

The assay of Benazepril and Amlodipine was performed with tablets and the % assay was found to be 100.08 and 100.66 which shows that the method is useful for routine analysis.

The linearity of Benazepril and Amlodipine was found to be linear with a correlation coefficient of 0.999 and 0.999, which shows that the method is capable of producing good sensitivity.

The acceptance criteria of precision is RSD should be not more than 2.0% and the method show precision 0.8 and 0.9 for Benazepril and Amlodipine which shows that the method is precise.

The acceptance criteria of intermediate precision is RSD should be not more than 2.0% and the method show

precision 0.5 and 0.7 for Benazepril and Amlodipine which shows that the method is repeatable when performed in different days also.

The accuracy limit is the percentage recovery should be in the range of 97.0% - 103.0%. The total recovery was found to be 100.42% and 100.64% for Benazepril and Amlodipine. The validation of developed method shows that the accuracy is well within the limit, which shows that the method is capable of showing good accuracy and reproducibility.

The acceptance criteria for LOD and LOQ are 3.00 and 10.00. The LOD and LOQ for Benazepril was found to be 3.04 and 10.02 and LOD and LOQ for Amlodipine was found to be 3.00 and 10.06.

The robustness limit for mobile phase variation and flow rate variation are well within the limit, the % degradation results are in limits. Which shows that the method is having good system suitability and precision under given set of conditions.

**BIBLIOGRAPHY**

1. Gurudeep Chatwal and Sham Anand, Instrumental methods of chemical analysis. Himalaya publishers, 7<sup>th</sup> edition, 1992; 2.624-2.639.
2. Skoog et al., Principles of Instrumental Analysis. Barkhanath publishers, 8<sup>th</sup> edition, 973-995.
3. Hobart. H. Willard et al., Instrumental methods of analysis, CBS Publications and Distributors, New Delhi, 1<sup>st</sup> edition, 1986, 529-563.
4. Sethi P.D., Quantitative analysis of Drugs & Pharmaceuticals. CBS publishers and distributors, New Delhi, 3<sup>rd</sup> edition, 2001; 1-120.
5. Janeyulu.Y & Marayyah, Quality Assurance & Quality Management in Pharmaceutical Industry. Pharma book publishers,Hyd, edition, 2005; 78-108.
6. Vogel's Text book of quantitative chemical analysis. Published by Dorling Kindersley pvt.ltd. 6<sup>th</sup> edition, 289-304.
7. Lloyd R. Snyder et al., practical HPLC method development. John Wiley & sons publishers, 2<sup>nd</sup> edition, 350-400.
8. Knevel A.M. & Digengl F.E, Jenkins Quantitative Pharmaceutical Chemistry, Mc Graw Hill Book Co.
9. Daniel W.Armstrong, Bonded Phase material for Chromatographic separations, 1985, U.S. Patent 4539399.
10. Sastry, C.S.P., Singh, N.R., Reddy, Methods of Analysis, 1986 edition, 316.